

Ш

ш

I

ഗ

4

CA

Z I

Ш

TECHNICAL DATA SHEET

ETHYLEN OXIDE TAPE

PRODUCT CODE : 4AEOT1950 / 4AEOT2550

RELATED STANDARDS: 4A Medical Ethylene Oxide Tape is classified as Type 1 according to ISO

11140-1 standard.

INTENDED USE : 4A Medical Ethylene Oxide Tapes are used to close wraps of the medical

devices to be sterilized and provide visual confirmation that the package

has been processed in a Ethylen Oxide sterilization.

OPERATION CONDITIONS: 4A Medical Ethylene Oxide Tapes are designed to operate on %60±:%10

humidity, 37°C±0,5°C, 90 min or 54°C±0,5°C 30 min at 600±30 mgratingg/L Ethylene Oxide concentration. Tape is latex free.

Before Sterilization

After Sterilization



The Ethylene Oxide Tape is not classified as dangerous.

PRODUCT : 19 mm x 50 m and 25 mm x 50 m Ethylene Oxide indicator printed on **SPECIFICATIONS** crepe paper tape.

TAPE	SPECIFICATION	UNIT	STANDARD	TYPICAL VALUE
	WIDTH	mm	-	19±1
	LENGTH	m	-	50± % 1
	THICKNESS	micron	1	110±10
	ADHESIVE	g/m^2	-	28±3
	INSTANT TACK	N/25 mm	ASTM D3121	350±50
	PEEL ADHESION	g/25 mm	ASTM D3121	450±80
	INDICATOR CONTROL	-	ISO 11140-1	Must return the specified color

STORAGE AND DISPOSAL After wrapping process tape can be used to prevent crepe and SMS paper

opening. Do not use damaged products. Keep products at room temperature 25°C. Shelf life is 4 years after production date.

PACKAGING: 19 mm x 50 m, 50 pieces in carton box / 25 mm x 50m, 40 pieces in

carton box